


Alkyne-PEG4-maleimide

<http://hk.lumiprobe.com/p/alkyne-peg4-maleimide>

When proteins or other sensitive biomolecules need to be conjugated, polyethyleneglycol (PEG) linkers are a great choice. They provide hydrophilicity and maintain separation between conjugated molecules. This linker based on PEG4, has length of 14Å (1.4 nm). Heterobifunctional linkers carrying different reactive groups can be used to prepare cross-conjugates between different molecules.

Maleimide functional group is especially useful, because they can label proteins site specifically. Most proteins have only limited number of cysteine residues that selectively react with maleimides. Alkyne group can be then conjugated with azides in a selective copper-catalyzed process, CuAAC (copper-catalyzed alkyne-azide cycloaddition).



外观:

质谱 M+ 增量: 311.3

分子量: 311.33

CAS 编号: 1262681-30-0

分子式: C₁₅H₂₁NO₆

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。