

## **LumiMount® Fluorescence Mounting Medium**

<http://hk.lumiprobe.com/p/lumimount-mounting-medium>

LumiMount® is an aqueous-based mounting medium for coverslipping cells and tissues stained with various fluorescent dyes.

LumiMount is compatible with most fluorescent markers used in microscopy applications and can reduce fluorochrome quenching during imaging due to anti-fade agent content. The refractive index of liquid mounting medium is 1.38, but it can increase to 1.46 as the water evaporates at the edges of the coverslip.

LumiMount provides a semi-permanent seal for prolonged slide storage at 2°C to 8°C. Since LumiMount is water-soluble, the coverslip may be removed by submerging the slide in a PBS solution until the coverslip is loosened. For better sealing of slides, it is recommended to seal the edges of the coverslip with waterproof nail polish.

---

外观:

溶解度:

质量控制:

储存条件: 儲存: 收到後室溫避光保存24個月。不要冷藏。避免暴露在光線下。

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。