

LumiMount® Plus Fluorescence Mounting Medium

<http://hk.lumiprobe.com/p/lumimount-plus-mounting-medium>

LumiMount® Plus is a ready-to-use aqueous-based mounting medium for coverslipping and imaging fluorescently labeled cell or tissue samples.

LumiMount Plus has a viscosity that provides easy application and uniform spreading over the tissue section, and it is not prone to bubbles.

Due to its unique formulation, LumiMount Plus reduces the rapid photobleaching of fluorescent proteins and dyes and prolongs possible imaging time. The refractive index of liquid LumiMount Plus is 1.43, but this property increases up to 1.52 after water evaporation and solidification of mounting.

LumiMount Plus provides a semi-permanent seal for prolonged slide storage at 2° to 8°C. Since LumiMount Plus is water-soluble, the coverslip may be removed by submerging the slide in a PBS solution until the coverslip is loosened.

LumiMount Plus is not compatible with intercalating dyes such as Hoechst, DAPI etc.

外观:

溶解度:

质量控制:

储存条件: 儲存: 收到後室溫避光保存24個月。不要冷藏。避免暴露在光線下。

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。