



**Lumiprobe Corporation**  
115 Airport Dr Suite 160  
Westminster, Maryland 21157  
美國  
手機: +1 888 973 6353  
傳真: +1 888 973 6354  
電子郵件: [order@lumiprobe.com](mailto:order@lumiprobe.com)

---

## Proteinase K Dilution Buffer

<http://hk.lumiprobe.com/p/proteinase-k-dilution-buffer>

Buffer solution for dissolving lyophilized proteinase K, maintaining pH 7.5. The buffer contains calcium chloride, which increases the thermal stability of proteinase K.

The recommended concentration of proteinase K stock solution is 20 mg/mL. After dilution, the proteinase K solution must be stored at a temperature of -25 to -18 °C.

This buffer is included in kits for DNA isolation [Spin Kit for Blood and Buccal Swab](#) and [Spin Kit for Any Sample](#).

---

外观:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。