

QuDye® dsDNA HS Reagent, 200×

<http://hk.lumiprobe.com/p/qureader-dsdna-hs-reagent>

QuDye® dsDNA HS Reagent 200×, quantification reagent is extremely selective for dsDNA over ssDNA and RNA. To prepare the working solution the concentrated QuDye dsDNA HS Reagent has to be diluted 1:200.

QuDye dsDNA HS Reagent usage is a high sensitive method over other traditional of DNA quantitation and is intended for use with all models fluorometers.

The 200× formulation is also a component of the nucleic acid quantification [QuDye dsDNA HS Assay Kit](#) for dsDNA samples measurements.

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。