

QuDye dsDNA HS Reagent, 200×

<http://hk.lumiprobe.com/p/qureader-dsdna-hs-reagent>

QuDye dsDNA HS Reagent 200×, quantification reagent is extremely selective for dsDNA over ssDNA and RNA. To prepare the working solution the concentrated QuDye dsDNA HS Reagent has to be diluted 1:200.

QuDye dsDNA HS Reagent usage is a high sensitive method over other traditional of DNA quantitation and is intended for use with all models fluorometers.

The 200× formulation is also a component of the nucleic acid quantification [QuDye dsDNA HS Assay Kit](#) for dsDNA samples measurements.

法律声 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經
明： 明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。