


Aldosterone-d4

<http://hk.lumiprobe.com/p/aldosterone-d4>

Aldosterone is a steroid hormone synthesized in the adrenal cortex. It plays a key role in regulating water-salt balance and blood pressure through its involvement in the renin-angiotensin-aldosterone system (RAAS). Aldosterone is used as a clinical marker in the diagnosis of primary aldosteronism (PA), as well as for assessing RAAS activity in renal disease, heart failure, and other conditions associated with water-electrolyte imbalance.

The product is intended for research purposes only.



外观: 白色粉末

分子量: 364.47

CAS 编号: 2368815-22-7

分子式: $C_{21}H_{24}D_4O_5$

溶解度:

质量控制: NMR 1H 和HPLC-MS (95+ & nbsp;% , d: 98+ & nbsp;%)

储存条件: 收到後 -20°C 避光保存 24 個月。 運輸: 室溫最多可保存3週。 乾燥。

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。