


Alkyne dT phosphoramidite

<http://hk.lumiprobe.com/p/alkyne-amidite-dt>

This product will be discontinued soon in favor of [Alkyne dT phosphoramidite](#). We may still have some material in stock (please see availability above).

Alkyne dT phosphoramidite is used for the synthesis of oligonucleotides containing an internal alkyne modification attached to the uridine 5-position. The alkyne group can be subsequently modified by various azides using copper-catalyzed Click chemistry.



外观:

分子量: 935.01

分子式: $C_{50}H_{59}N_6O_{10}P$

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。

偶联条件:

切割条件:

解保护条件: