


Alkyne CPG 1000A

<http://hk.lumiprobe.com/p/alkyne-cpg-1000>

Controlled pore glass solid support for synthesizing 3'-alkyne modified oligonucleotides. The support allows the preparation of oligos of up to 100 base pairs. The modification is stable in standard deblocking conditions and can tolerate ammonia, AMA, and other basic reagents.



外观:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。

孔径大小, 埃: 1000

典型載荷, $\mu\text{mol/g}$: 40–70

偶聯條件:

切割條件: