


---

## Alkyne-PEG3-amine

<http://hk.lumiprobe.com/p/alkyne-peg3-amine-932741-19-0>

Alkyne-PEG3-amine is a heterobifunctional molecule based on triethyleneglycol (PEG3) scaffold. It contains amine and alkyne functional groups. The former can be conjugated with various electrophiles. The latter is suitable for copper-catalyzed click chemistry and Sonogashira coupling.

PEG linkers are hydrophilic linkers of choice for conjugating biomolecules.



---

外观:

分子量: 187.24

CAS 编号: 932741-19-0

分子式: C<sub>9</sub>H<sub>17</sub>NO<sub>3</sub>

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。