


## Alkyne-PEG3-iodide

<http://hk.lumiprobe.com/p/alkyne-peg3-iodide>

Alkyne-PEG3-iodide (or Iodo-PEG3-alkyne) is a hydrophilic bifunctional PEG derivative for the alkylation of nucleophilic groups and copper-catalyzed click chemistry. As compared with chloride and bromide, iodide is the halogen group that is the most active in alkylation reactions.



外观:

分子量: 298.12

CAS 编号: 1234387-35-9

分子式: C<sub>9</sub>H<sub>15</sub>IO<sub>3</sub>

IUPAC 名称: 3-[2-[2-(2-iodoethoxy)ethoxy]ethoxy]-1-propyne

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。