


Alkyne-PEG4-amine

<http://hk.lumiprobe.com/p/alkyne-peg4-amine-1013921-36-2>

Alkyne-PEG4-amine is a bifunctional oligoethylene glycol linker containing both an amine group for the easy attachment to various electrophilic reagents (NHS esters, epoxides, etc.), and an alkyne group for click chemistry. This allows a smooth conjugation of different biomolecules.

Reasonably long, highly hydrophilic PEG4 linker ensures smooth conjugation and separation of conjugated fragments in aqueous solutions.



外观:

分子量: 231.29

CAS 编号: 1013921-36-2

分子式: C₁₁H₂₁NO₄

IUPAC 名称: 3,6,9,12-Tetraoxapentadec-14-yn-1-amine

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。