


Androstenedione-d7

<http://hk.lumiprobe.com/p/androstenedione-d7>

Androstenedione is a steroid hormone produced in testes in the ovary and in the adrenal cortex. Androstenedione is a precursor of testosterone and estrone. Androstenedione originates either from the conversion of dehydroepiandrosterone or from 17-hydroxyprogesterone. Androstenedione (ASD) is a biomarker used in the diagnostic workup of hyperandrogenism, congenital adrenal hyperplasia, premature adrenarche, and polycystic ovary syndrome (PCOS).



外观:

分子量: 293.45

CAS 编号: 67034-85-9

分子式: $C_{19}H_{19}D_7O_2$

质量控制: NMR 1H and HPLC-MS (95+ %, D: 98+ %)

储存条件: 收到後在 $-20^{\circ}C$ 避光條件下保存 24 個月。運輸: 室溫最多可保存3週。乾燥。

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。