

3-Azidopropanol

<http://hk.lumiprobe.com/p/azidopropanol-72320-38-8>

Azidopropanol is a bifunctional synthone carrying azide and hydroxyl groups. It can be used for the click chemistry or Staudinger ligation involving its azide group.

Hydroxyl group can be functionalized in a non-aqueous environment by means of reaction with CDI and subsequent coupling with amines producing carbamates.

外观:

分子量: 101.11

CAS 编号: 72320-38-8

分子式: $C_3H_7N_3O$

IUPAC 名称: 3-azidopropanol-1

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。