


## BDP® 493/503 lipid stain

<http://hk.lumiprobe.com/p/bdp-493-503-lipid-stain>

BDP 493/503 is a green fluorescent borondipyrromethene dye for lipid labeling. It is used for cell membrane visualization in microscopy, as well as for staining and tracing of hydrophobic media.



外观:

分子量: 262.11

CAS 编号: 121207-31-6

分子式: C<sub>14</sub>H<sub>17</sub>N<sub>2</sub>BF<sub>2</sub>

IUPAC 名称: [2-[1-(3,5-dimethyl-2H-pyrrol-2-ylidene)ethyl]-3,5-dimethyl-1H-pyrrolato-N1,N2]difluoroboron

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。