


## Biotin-PEG4-carboxylic acid

<http://hk.lumiprobe.com/p/biotin-peg4-carboxylic-acid>

Biotin-PEG4-carboxylic acid is a bifunctional tetraethyleneglycol derivative bearing a carboxy group and biotin moiety. The carboxy group can be activated with peptide coupling reagents like PyBOP or carbodiimides like EDC to form a stable amide linkage with amines. Biotin-labeled biomolecules can be bound to avidin or streptavidin for further purification and detection.

The structure of this reagent features a long hydrophilic PEG4 linker that separates biotin residue from the target molecule for efficient binding with streptavidin. The linker also enhances aqueous solubility to facilitate conjugation.



外观:

分子 491.61

量:

CAS 721431-18-1

编号:

分子 C<sub>21</sub>H<sub>37</sub>N<sub>3</sub>O<sub>8</sub>S

式:

IUPAC 3-[2-[2-[2-[5-[(3aS,4S,6aR)-2-oxo-1,3,3a,4,6,6a-hexahydrothieno[3,4-d]imidazol-4-yl]pentanoylamino]ethoxy]ethoxy]ethoxy]propanoic  
名称: acid

溶解度:

质量控制:

储存条件:

法律 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何  
聲明: 用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。