


(C26) Carnitine-d3

<http://hk.lumiprobe.com/p/c26-carnitine-d3>

(C26) Carnitine-d3 is a deuterium labeled (C26) Carnitine and intended for use as an internal standard in GC- or LC-MS.

(C26) D,L-Carnitine could be used as a diagnostic marker in dried blood spots from newborns and patients with X-linked adrenoleukodystrophy. (C26) D,L-Carnitine could also be a marker for Zellweger spectrum disorders (ZSD) which are a group of genetic metabolic disorders caused by a defect in peroxisome biogenesis.



外观:

分子量: 579.37

分子式: $C_{33}H_{63}D_3ClNO_4$

质量控制: NMR 1H 和 HPLC-MS (95+%, D: 98+%)

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。