


## C26:0-lysophosphatidylcholine (LPC)

<http://hk.lumiprobe.com/p/c26-lpc-1213783-80-2>

Lysopc(26:0), also known as LysoPC a C26:0, is classified as a member of the 1-acyl-sn-glycero-3-phosphocholines. 1-acyl-sn-glycero-3-phosphocholines are glycerophosphocholines in which the glycerol is esterified with a fatty acid at O-1 position, and linked at position 3 to a phosphocholine. Lysopc(26:0) can be used as internal standard LC-MS studies.

Lysophosphatidylcholines (LPCs) specifically bind to G protein-coupled receptors GPR119, GPR40, GPR55 and GPR4. Binding of LPCs to GPR119, GPR40 and GPR55 induces intracellular calcium mobilization and leads to increased glucose-stimulated insulin secretion in different cell systems. In blood or plasma LPCs are bound mainly to albumin and to a lesser extent to lipoproteins.

Inflammation, cell damage and other pathophysiological conditions can profoundly alter the ratio of free to albumin bound LPC through increased production of LPC or decreased plasma levels of albumin. In particular, lower levels of albumin (hypoalbuminemia) lead to lower levels of LPC in the blood. Lysopc(26:0) is a biomarker for X-linked Adrenoleukodystrophy (X-ALD). Lysopc(26:0) could be used as auxiliary marker for diagnosis of peroxisomal  $\beta$ -oxidation disorders.



外观: 白色粉末

分子量: 635.90

CAS 编号: 1213783-80-2

分子式:  $C_{34}H_{70}NO_7P$

IUPAC 名称: [(2R)-3-hexacosanoyloxy-2-hydroxypropyl] 2-(trimethylazaniumyl)ethyl phosphate

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。