


DOTA-PEG4-amine hydrochloride

<http://hk.lumiprobe.com/p/dota-peg4-amine>

DOTA-PEG4-amine contains a terminal amino group, DOTA (1,4,7,10-Tetraazacyclododecane-1,4,7,10-tetraacetic acid) moiety, and a linear PEG chain as a linker. It is widely used for bioconjugation in the imaging diagnostics field due to the DOTA possessing chelating properties, and a highly reactive amine group.

The hydrophilic PEG chain increases the aqueous solubility. PEG linker allows to incorporate the chelating agent by conjugation of DOTA-PEG4-amine with carboxylic acids and activated esters using even sterically hindered moieties in a biomolecule.



外观:

分子量: 760.97

分子式: $C_{24}H_{51}N_6Cl_5O_{10}$

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。