


## DOTA-PEG4-azide hydrochloride

<http://hk.lumiprobe.com/p/dota-peg4-azide-hydrochloride>

DOTA (dodecane tetraacetic acid, also known as tetraxetan) is a complexing agent, especially for lanthanide ions. DOTA-PEG4-azide is a bifunctional linker containing chelating and azide groups for conjugation with alkyne-containing biomolecules. The structure of the reagent features a long hydrophilic PEG4 linker that separates the DOTA residue from the target molecule. The linker also enhances aqueous solubility to facilitate conjugation.

DOTA-PEG4-azide can be used to label radiotherapeutic agents or imaging probes and detect labeled tumors with PET, SPECT, and CT methods.



外观:

分子量: 750.51

分子式: C<sub>24</sub>H<sub>48</sub>Cl<sub>4</sub>N<sub>8</sub>O<sub>10</sub>

溶解度:

质量控制:

储存条件:

法律声明: