


## dsGreen®, solid form

<http://hk.lumiprobe.com/p/dsgreen-i-solid-form>

dsGreen® is a highly sensitive and selective stain for double-stranded DNA (dsDNA). dsGreen solutions are used for DNA binding and [gel staining](#). This product is a solid form of dsGreen that can be used to prepare solutions according to research protocols. Lumiprobe also manufactures dsGreen in the form of [a solution](#) specifically optimized for real-time PCR.



外观:

分子量: 509.73

CAS 编号: 163795-75-3

分子式:  $C_{32}H_{37}N_4S$

溶解度:

质量控制:

储存条件:

法律声明:

本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。

激发/吸收极大值, 纳米: 454

发射极大值, 纳米: 524