



## FAM tetrazine, 5-isomer

<http://hk.lumiprobe.com/p/fam-tetrazine-5>

Tetrazine (methyltetrazine) derivative of fluorescein, pure 5-isomer.

Tetrazines are reactive heterodiene components for [4+2] inverse electron demand Diels-Alder reaction (ied-DA) with strained cycloalkenes (trans-cyclooctenes, cyclopropenes). This reaction is also called TCO (tetrazine-cyclooctene) ligation. The reaction is among the fastest bioconjugation reactions.



外观:

质谱 M+ 增量: 529.12

分子量: 559.53

CAS 编号: 1619222-85-3

分子式: C<sub>31</sub>H<sub>21</sub>N<sub>5</sub>O<sub>6</sub>

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 492

ε, 摩尔吸光系数, cm<sup>-1</sup>: 74000

发射极大值, 纳米: 517

荧光量子产率: 0.93

CF<sub>260</sub>: 0.22

CF<sub>280</sub>: 0.17