

## **QuReader® 8 Green**

<http://hk.lumiprobe.com/p/fluorometer-qureader-8-well-green>

The QuReader® 8 Green fluorometer is a desk fluorometer with a touch screen designed for dsDNA, ssDNA, RNA, and protein quantitation. The QuReader® fluorometer enables the fluorescence measuring for up to eight samples using green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acid and protein concentrations in the lower pg/uL range.

---

法律聲明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。