

QuReader® 8 Green

<http://hk.lumiprobe.com/p/fluorometer-qureader-8-well-green>

The QuReader® 8 Green fluorometer is a desk fluorometer with a touch screen designed for dsDNA, ssDNA, RNA, and protein quantitation. The QuReader® fluorometer enables the fluorescence measuring for up to eight samples using green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acid and protein concentrations in the lower pg/uL range.

法律声 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經
明： 明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。