

Iodo-PEG3-carboxylic acid

<http://hk.lumiprobe.com/p/i-peg3-carboxylic-acid>

A molecule containing carboxylic acid (COOH) group and alkyl iodide functional group.

Alkyl iodides are powerful alkylators. They react with a variety of nucleophiles, very active against thiol groups, amines, alcoholates, and phenolates.

This molecule allows performing a variety of conjugations via a flexible and hydrophilic PEG3 linker.



外观:

分子量: 332.13

分子式: $C_9H_{17}IO_5$

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。