


Iodo-PEG3-carboxylic acid

<http://hk.lumiprobe.com/p/i-peg3-carboxylic-acid>

A molecule containing carboxylic acid (COOH) group and alkyl iodide functional group.

Alkyl iodides are powerful alkylators. They react with a variety of nucleophiles, very active against thiol groups, amines, alcoholates, and phenolates.

This molecule allows performing a variety of conjugations via a flexible and hydrophilic PEG3 linker.



外观:

分子量: 332.13

分子式: C₉H₁₇IO₅

溶解度:

质量控制:

储存条件: