

## **Magnetic Beads (100 mg/mL)**

<http://hk.lumiprobe.com/p/magnetic-beads>

Magnetic silica-based particles represent the modern standard for rapid and efficient DNA and RNA isolation. This technology delivers preparations of the highest purity, free from inhibitors and with minimal sample loss. Offering seamless compatibility with automated platforms, it is equally suited for manual protocols and high-throughput research. You will obtain pure, ready-to-use nucleic acids, ideal for even the most demanding downstream applications.

---

法律聲明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。