



## LumiTracker® Mito Red CMXRos

<http://hk.lumiprobe.com/p/mitotracker-red-cmxros>

LumiTracker Mito Red CMXRos is a cationic, red-fluorescent dye for staining active mitochondria in live cells. CMXRos is a hydrophobic compound that passively diffuses across the plasma membrane and selectively accumulates in active mitochondria based on their membrane potential. LumiTracker Mito Red CMXRos is compatible with aldehyde fixation. This dye is useful to indicate cell health as well as for mitochondria localization.

Ros in the name of this product stands for rosamine and has nothing to do with the detection of reactive oxygen species (ROS).



外观:

分子量: 531.53

CAS 编号: 167095-09-2

分子式: C<sub>32</sub>H<sub>32</sub>Cl<sub>2</sub>N<sub>2</sub>O

溶解度:

质量控制:

储存条件:

法律声明:

本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。

激发/吸收极大值, 纳米: 581

ε, 摩尔吸光系数, cm<sup>-1</sup>: 101000

发射极大值, 纳米: 600

荧光量子产率: 0.91