

---

## Random(dN)6 primer

<http://hk.lumiprobe.com/p/n6-stock-oligonucleotide>

A random hexamer or random-sequence oligonucleotide can be used as a universal primer for cDNA synthesis by reverse transcription and reverse transcription PCR (RT-PCR) technique. One of the commonly used oligonucleotides.

## Dilution

Original quantity, nmol	Final concentration, uM	Volume of water needed to dilute a solution to a desired concentration, mL
	100	0.5
50	50	1
	20	2.5

---

外观:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。