


---

## PEG3-dicarboxylic acid

<http://hk.lumiprobe.com/p/peg3-dicarboxylic-acid-31127-85-2>

PEG3-dicarboxylic acid is a homobifunctional linker with two carboxylic acid groups flanking the PEG3 molecule.

The linker is highly hydrophilic and useful for conjugating and immobilizing proteins, peptides, and other biomolecules. The carboxy groups can be activated by carbodiimides.



---

分子量: 294.3

CAS 编号: 31127-85-2

分子式: C<sub>12</sub>H<sub>22</sub>O<sub>8</sub>

IUPAC 名称: 3-[2-[2-(2-(2-carboxyethoxy)ethoxy)ethoxy]ethoxy]propanoic acid, 4,7,10,13-tetraoxahexadecanedioic acid

溶解度:

质量控制:

储存条件: