


PEG3-dicarboxylic acid

<http://hk.lumiprobe.com/p/peg3-dicarboxylic-acid-31127-85-2>

PEG3-dicarboxylic acid is a homobifunctional linker with two carboxylic acid groups flanking the PEG3 molecule.

The linker is highly hydrophilic and useful for conjugating and immobilizing proteins, peptides, and other biomolecules. The carboxy groups can be activated by carbodiimides.



分子量: 294.3

CAS 编号: 31127-85-2

分子式: C₁₂H₂₂O₈

IUPAC 名称: 3-[2-[2-[2-(2-carboxyethoxy)ethoxy]ethoxy]ethoxy]propanoic acid, 4,7,10,13-tetraoxahexadecanedioic acid

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。