


PEG3-tetrazine

<http://hk.lumiprobe.com/p/peg3-tetrazine>

PEG3-tetrazine is a monoreactive derivative that allows a facile attachment of hydrophilic triethyleneglycol fragment (pegylation) by inverse electron demand Diels-Alder reaction (IEDDA) with strained olefins, such as trans-cyclooctenes, and some cyclooctynes.



外观:

分子量: 334.37

分子式: C₁₆H₂₂N₄O₄

IUPAC 名称: 2-(2-{2-[2-(6-Phenyl-1,2,4,5-tetraazin-3-yl)ethoxy]ethoxy}ethoxy)ethanol

溶解度:

质量控制:

储存条件:

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。