


## Phosphate CPG 1000 solid support

<http://hk.lumiprobe.com/p/phosphate-cpg-1000>

Phosphate modified oligos are used for gene synthesis by ligation, and other applications. Phosphate CPG 1000 is a modifier solid support for the synthesis of oligonucleotides bearing phosphate at 3 prime end.

CPG 1000 is suitable for the synthesis of long oligos up to 120 mer, the use of long oligos allows for a more efficient gene assembly. The CPG is compatible with standard synthesis and deblocking procedures, it ensures smooth cleavage and complete phosphorylation.



外观:

质量控制:

储存条件:

孔径大小, 埃: 1000

典型载荷,  $\mu\text{mol/g}$ : 30–50

偶联条件:

切割条件:

解保护条件: