


Phosphate CPG 1000 solid support

<http://hk.lumiprobe.com/p/phosphate-cpg-1000>

Phosphate modified oligos are used for gene synthesis by ligation, and other applications. Phosphate CPG 1000 is a modifier solid support for the synthesis of oligonucleotides bearing phosphate at 3'-end.

CPG 1000 is suitable for the synthesis of long oligos up to 100 base pairs, the use of long oligos allows for a more efficient gene assembly. The CPG is compatible with standard synthesis and deblocking procedures, it ensures smooth cleavage and complete phosphorylation.



外观:

质量控制:

储存条件:

法律声明:

本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試，且未經明示或暗示授權用於其他任何用途，包括但不限於體外診斷、人類或動物用途，以及商業用途。

孔径大小, 埃: 1000

典型载荷, $\mu\text{mol/g}$: 30–50

偶联条件:

解保护条件: