



## Pyrenebutyric acid

<http://hk.lumiprobe.com/p/pyrenebutyric-carboxylic-acid>

Pyrenylbutyric acid is a derivative of pyrene hydrocarbon with a free carboxylic acid function. The reagent is useful as a control for experiments with other reactive pyrene derivatives such as [pyrene NHS ester](#) and [pyrene azide](#). Carboxylic acid function can be activated by carbodiimides and other activating reagents.



外观:

分子量: 288.34

CAS 编号: 3443-45-6

分子式:  $C_{20}H_{16}O_2$

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 343; 326; 313; 276; 265; 242; 234

发射极大值, 纳米: 377; 397