
DryDrops® qPCR, Single Assay qPCR Lyophilized Beads

<http://hk.lumiprobe.com/p/qpcr-lyophilized-beads>

Single Assay qPCR Lyophilized Beads are a pre-formulated and pre-dispensed reaction mixture for carrying out polymerase chain reaction (PCR). Each lyophilized bead contains all the necessary components for PCR and qPCR with a volume of 25 µL:

- hot-start Taq polymerase;
- mixture of deoxynucleoside triphosphates (dNTPs);
- optimized buffer for PCR (contains Mg²⁺ with a concentration of 3 mM in 1× reaction mixture);
- protectors for lyophilization.

The composition of the reaction mixture is optimized to obtain ideal results in terms of processivity and specificity of amplification. Lyophilized beads provide excellent reproducibility between reactions by minimizing pipetting steps and reducing the likelihood of pipetting errors and sample contamination. Each batch of the Single Assay qPCR Lyophilized Beads undergoes functional testing to ensure batch-to-batch reproducibility.

The beads are suitable for real-time PCR detection, including quantitative analysis (with fluorescent probes or intercalating dye, e.g., [Eva488](#)), as well as for DNA amplification with subsequent detection of results by electrophoresis. The beads are also suitable for routine cloning tasks and other applications that require further handling of the PCR product after amplification.

Single Assay qPCR Lyophilized Beads are compatible with any type of amplifier and can be used for PCR in amplifiers with classic thermal blocks for PCR tubes and amplifiers with cartridges.

The lyophilized format allows storing the PCR mixture for up to 12 months at temperatures up to 4 °C.

We offer Single Assay qPCR Lyophilized Beads in vials of 10 pcs. (**11534**), vials of 50 pcs. (**21534**), and PCR-strips of 8 pcs. (**31534**).

外观: 白色光滑球體

法律声明: 本產品僅供研究目的提供和銷售。本產品並未經過食品、藥品、醫療器械、化妝品等領域的安全性和效力測試, 且未經明示或暗示授權用於其他任何用途, 包括但不限於體外診斷、人類或動物用途, 以及商業用途。