



## sulfo-Cyanine7 bis-NHS ester

<http://hk.lumiprobe.com/p/sulfo-cy7-bis-nhs-ester>

sulfo-Cyanine7 bis-NHS ester is a water-soluble, bifunctional cross-linker containing a core of near-infrared, water-soluble sulfo-Cyanine7, and two amine-reactive NHS ester groups.



外观:

分子量: 1041.23

分子式: C<sub>50</sub>H<sub>57</sub>N<sub>4</sub>KO<sub>4</sub>S<sub>2</sub>

溶解度:

质量控制:

储存条件:

法律声明:

激发/吸收极大值, 纳米: 750

ε, 摩尔吸光系数, cm<sup>-1</sup>: 240600

发射极大值, 纳米: 773

荧光量子产率: 0.24

CF<sub>260</sub>: 0.04

CF<sub>280</sub>: 0.04