

## **TechnoDye 502/512**

<http://hk.lumiprobe.com/p/technodye-502-512>

TechnoDye 505/512 is a hydrophobic technical dye with a bright fluorophore with an excitation spectrum in the green range. Photophysical properties, such as strong absorption in the visible region, high fluorescence quantum efficiency, and high photostability, make it suitable for different applications.

---

外观:

分子量: 248.08

分子式:  $C_{13}H_{15}BF_2N_2$

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 505

发射极大值, 纳米: 512